Instruments

MAPN

Mikroskopisch assistierte perkutane Nukleotomie
Microscopically Assisted Percutaneous Nucleotomy
Nucleotomía percutánea microscópicamente asistida
Nucléotomie percutanée assistée au microscope
Nucleotomia percutanea assistita al microscopio

<https://ra.medicon.de>

Not all products illustrated in this documentation are available for sale in all countries. The products shown in this catalogue comply with the Medical Device Directive 93/42/EEC (MDD).

Just scan the QR-Code
and you already
reach the latest
information!

medicon
Instrumente
the art of surgery

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With an incidence rate of 87 per 100 000 inhabitants per year in Germany [1], lumbar intervertebral disk surgery is one of the most common procedures in the fields of Neurosurgery and Orthopaedics. This high frequency rate requires procedures that, on the one hand, afford extensive, qualitative safety during surgery and, on the other hand, keep access-related trauma to a minimum, thus promoting rapid patient rehabilitation. Healthcare costs can be reduced by shortening the length of time spent in hospital and through the availability of an alternative, minimally invasive intervertebral disk surgical technique performed on an outpatient basis.

The previous open intervertebral disk procedure was replaced in the 1970s by the microsurgical technique developed by Yasargil [2] and Caspar [3], and which nowadays is considered as the "Gold standard". The advantage of this method compared with the "open" technique is the smaller access required and the excellent visualisation of the surgical field using a surgical microscope. Another important step towards reducing access trauma was achieved with the development of the endoscopically supported method by Destandeu [4] and Foley [5] for the surgical correction of lumbar herniated vertebral disks. The disadvantage of this method mainly lies in the two-dimensional depiction of the surgical field and the troublesome contamination of the optics due to bleeding.

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[2] **YASARGIL** Microsurgical operation of the herniated lumbar disc. In: WÜLLENWEBER, BROCK, HAMER (eds.) Advances in Neurosurgery, Springer-Verlag, Berlin 1977

[3] **CASPAR** A new surgical procedure for lumbar disc herniation causing less tissue damage through a microsurgical approach. In: WÜLLENWEBER, BROCK, HAMER (eds.) Advances in Neurosurgery. Springer-Verlag, Berlin 1977

[4] **DESTANDEAU** A special device for endoscopic surgery of lumbar disc herniation. Neurol Res 21 (1999) 39-41

[5] **FOLEY, SMITH** Microendoscopic discectomy. Tech Neursurg 3 (1997) 301-307

Mikroskopisch assistierte perkutane Nukleotomie
Microscopically assisted percutaneous nucleotomy
Nucleotomía percutánea microscópicamente asistida
Nucléotomie percutanée assistée au microscope
Nucleotomia percutanea assistita con microscopio

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- ▶ This technique, referred to as microscopically assisted percutaneous nucleotomy (MAPN) allows optimal, three-dimensional visualisation of the surgical field through the surgical microscope.
- ▶ Access to the spinal canal is provided by a working channel, introduced through transmuscular dilatation.
- ▶ Without sacrifices in quality the same discal pathologies can be treated as with the standard microsurgical technique [6, 7].
- ▶ Neither wound drainage nor post-operative orthosis is required, thus enhancing patient comfort.
- ▶ Moreover the same technique allows a less invasive variant of access to all other levels of the posterior spinal column [7].

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MICROSCOPICALLY ASSISTED PERCUTANEOUS NUCLEOTOMY

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The patient leans on his/her knees and elbows, thus ensuring that the lumbar spine is not curved. A skin incision 15 mm long and approximately 2 cm from the midline is made at the level of the marker needle on the affected side. The level is marked with the needle on the uninvolved side, monitoring the lateral C-arm image. The needle should be directed towards the lesion (e.g. cranially or caudally extruded sequestrum). [1-3]

Insertion of the first dilatator and confirmation by lateral C-arm image. This prevents intervention at the wrong level. The target point is the Lig. flavum in the interlaminar window. Medial convergence is approximately 10°. The lower edge of the upper Lamina medial from the joint mass is a good anatomical landmark that is palpable with the dilatator. [4-5]

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

Insertion of the second dilatator with image intensifier control and removal of the first. Once again, the lower edge of the upper Lamina should be palpated as a safety precaution. The external screw allows the working channel to be "gently screwed in" through the paraspinal musculature above the second dilatator. Image intensifier control should then be repeated as a precautionary measure. [6-8]


The handle is placed on the working channel, thus facilitating straightforward positioning of the working channel in the direction of the affected site. Soft tissue is removed from the Lig. flavum beneath the microscope. Blunt opening of the Lig. flavum using a dissector. [9-10]

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
Resection of parts of the flavum to extend access into the spinal canal to reveal the dural tube and/or nerve root. Mobilisation of nerve root or dural sack. Protection of the neural structures using the nerve root retractor. Removal of the herniated disk. Removal or "unscrewing" of the working channel under observation in order to stem any bleeding from the paraspinal musculature. Examination beneath the surgical microscope shows how the muscle fibres overlie each other on retracting the working channel. There is no need to insert a suction tube. [11-14]

The wound is closed by one single suture of the fascia and one for the subcutis. Skin adaptation by a bandaid. The figure illustrates the cosmetic outcome 3 months after surgery. [15]

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© Copyright 12/2008, MEDICON eG, Tuttlingen