Instruments

MICRODISC XS

Miniaturisierte (XS) Mikrodiskektomie
Miniaturized (XS) Micro-Discectomy
Microdiscektomia miniaturizada (XS)
Mikrodiscectomie miniature (XS)
Mikrodiscectomie miniaturizzata (XS)

<https://ra.medicon.de>

Not all products illustrated in this documentation are available for sale in all countries. The products shown in this catalogue comply with the Medical Device Directive 93/42/EEC (MDD).

Just scan the QR-Code
and you already
reach the latest
information!

medicon
Instrumente
the art of surgery

Miniaturisierte (XS) Mikrodiskektomie
Miniaturized (XS) Micro-Discectomy
Microdiscectomía miniaturizada (XS)
Microdiscectomie miniature (XS)
Microdiscectomie miniaturizzata (XS)

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Depending on its location, a lumbar intraspinal disc herniation may be approached via an interlaminar or translaminar access. The paraspinal route is indicated for extraforaminal disc herniation. Herniated cervical, and more rarely, thoracic discs may also be approached via the posterior route, provided that at least 2/3 of the herniated volume is located lateral to the dural sac. In these cases, cervical dorsolateral foraminotomy and thoracic transfacetall approach are indicated.

At the end of the 1970s, independently of one another, Williams [1], Caspar [2], and Yasargil [3] introduced micro-surgery for the operative treatment of lumbar disc herniation. The advantages of minimal tissue trauma combined with early postoperative mobilization of the patients helped this method to become popular. As in conventional surgery, microdiscectomy requires dissection of the paravertebral muscles from the spinous process and hemilamina. Retraction due to scarring as well as denervation of the muscles, particularly rich in proprioceptors in the cervical region, may result in troublesome local syndromes.

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Endoscopic removal of disc fragments is performed through a cylindrical working channel (approx. 7mm Ø) introduced by muscle-splitting. The technique, developed by Destandeu [4] and Kambin [5] and refined by Foley [6] and Rütten [7], is very atraumatic to the muscles. However, its disadvantages are the two-dimensional view of the surgical field, the difficulty of using a burr and impairment of the optics by bleeding.

The Microdisc-XS technique combines the less traumatic "endoscopic" muscle-splitting approach with the proven "microsurgical" advantages in the surgical field. The spreading sleeve facilitates the introduction of three instruments, including a burr. Due to the reduced surgical traumatization the postoperative intake of analgesics has been significantly reduced [8].

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- [1] Williams (1978) Microlumbar discectomy. A conservative surgical approach to the virgin herniated disc. Spine, 3(2): 175-182
- [2] Caspar (1977) A new surgical procedure for lumbar disc herniation causing less tissue damage through a microsurgical approach. In: Advances in Neurosurgery; Wüllenweber, Brock, Hamer (eds); Vol 4, 74-77; Springer-Verlag, Berlin
- [3] Yasargil (1977) Microsurgical operation of the herniated lumbar disc. In: Advances in Neurosurgery; Wüllenweber, Brock, Hamer (eds); Vol 4: 81-88; Springer-Verlag, Berlin
- [4] Destandeu (1999) A special device for endoscopic surgery of lumbar disc herniation. Neurol Res 21: 39-41
- [5] Kambin, Zhou (1997) Arthroscopic discectomy of the lumbar spine. Clin Orthop Relat Res 337: 49-57
- [6] Foley, Smith, Rampersaud (1999) Microendoscopic approach to farlateral lumbar disc herniation. Neurosurg Focus 7(5): e5
- [7] Rütten, Komp, Godolias (2006) A new full-endoscopic technique for the interlaminar operation of lumbar disc herniations using 6-mm endoscopes: prospective 2-year results of 331 patients. Minim Invasive Neurosurg 49(2):80-7
- [8] Greiner-Perth R, Boehm H, Ezzati E, Hempel B, Leonhardi J, Franke J (2006) Comparison of soft tissue trauma in microsurgical nucleotomy to a new technique : Microscopically Assisted Percutaneous Nucleotomy. An MRI-Study. Eur Spine J 15 (Suppl.4):S475

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X-ray labelling: (lateral view) Paramedian insertion of the cannula contralateral to the disc herniation (DH) and projecting on the intervertebral space (IVS). If necessary, an additional marker can be placed at the level of the extruded disc fragment.

Skin incision: Under the microscope, approx. 18 mm long, 15 mm from the midline, centred on the IVS or DH. The surgeon's oculars are tilted 5°-10° towards the midline.

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Muscle splitting: With the index finger until contact is made with the lamino-facet junction (target area). Alternatively: sequential dilators. The target area is exposed with narrow Langenbeck retractors. A few insertions of the multifidus muscle are coagulated and divided.

Insertion of the spreading sleeve (SS): Replace the Langenbecks with Love nerve root retractors, and advance a SS of suitable length with alternating rotating movements until contact is made with bone. Alternatively: after removing the Langenbecks, advance the SS through a dilator as far as the target area.

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Interlaminar fenestration: Burr or punch out the lower border of the upper hemilamina and/or the medial part of the joint. Narrow burr-handpieces facilitate visual control. The yellow ligament is opened. The disc fragment is removed, and, if necessary, discectomy is performed by standard microsurgical technique. The nerve root retractor can be replaced by the flattened sucker in right or left version (3 sizes).

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The Trans-laminar Approach

Positioning: The lumbar hemilamina (HL) dips downwards obliquely in the caudocephalad direction. By tilting the operating table about 15° head upwards, the HL becomes horizontal. That facilitates drilling of the translaminar hole (TH).

X-ray labelling: (lateral view) Paramedian insertion of the cannula contralateral to the DH and projecting on the upper rim of the IVS. The second marker is placed at the level of the extruded disc fragment, usually between the lower border of the cranial pedicle and the upper rim of the IVS.

Skin incision: Under the microscope, approx. 18 mm long, 10 mm from the midline, centered on the DH. The surgeon's oculars are tilted 5°-10° towards the midline.

Fascia incision: Straight, a few millimetres longer than the skin incision, 2 retaining sutures.

Muscle dilatation: With the index finger until contact is made with the HL and its inferior border. Alternatively: sequential dilators. The HL is exposed with narrow Langenbeck retractors. Its lateral and inferior border as well its junction with the spinous process should be identified.




Insertion of the spreading sleeve (SS): Replace the Langenbecks with Love nerve root retractors and advance a SS of suitable length with alternating rotating movements until contact is made with the HL. Alternatively: after removing the Langenbecks, the SS is advanced through a dilator as far as the HL.

Fixing the SS: The SS docks the contralaterally fixed snake arm and is locked in the desired position. (1) Spread the SS. (2) Introduce a dissector into the SS and tampon with a moist unfolded sponge. (3) By means of lateral fluoroscopy adjust if necessary position and/or angle of the SS. The DH should project into the geometric centre of the oval opening of the SS.

Translaminar hole: The second skin marking, between the lower border of the cranial pedicle and the superior rim of the IVS, can be transferred to the HL using a small diamond burr. Centered on this, a TH about 10mm in diameter is drilled. Anatomical remarks: (a) The HLs become narrower in the caudocephalad direction, i.e. the TH turns from round to ovale; the IVS "ascends" increasingly under the HL and its superior rim may be inspected through the TH. (b) Keep a distance of at least 3 mm between the lateral border of the TH and the lateral rim of the HL in order to avoid a fracture of the pars interarticularis. (c) When drilling through the HL, there are 3 successive layers from the surface down: white (lamina externa), red (spongy bone), and white again (lamina interna). The latter should be drilled with a diamond burr as there is no protective yellow ligament beneath the lamina interna.

The TL should ideally be centered on the lateral border of the dural sac and the axilla of the exiting root. However, the root is often displaced cephalad by disc fragments. The DH is removed, and, if necessary, the IVS is cleared using standard microsurgical technique. Due to the small surgical field, the nerve root retractor can be replaced by the flattened sucker in right and left versions (3 sizes).

Wound closure: The SS is closed and withdrawn slowly. Any minor bleeding from muscle is coagulated. The access tunnel is filled with warm Ringer solution. 2-3 fascial sutures. 2 subcutaneous sutures. Intracutaneous suture or adhesive strips.

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




The Extra-foraminal Approach

Positioning: According to surgeon's preference (Wilson frame; knee-elbow position).

X-ray labelling: (lateral view) Paramedian insertion of the cannula contralateral to the DH and projecting on the IVS. AP-view: A rectangle is formed by 2 horizontal lines: (1) upper rim of the IVS (2) line projected on the middle third of the cranial transverse processes and by 2 vertical lines: (3) paramedian to the spinous process (4) lateral parapedicular line. The lateral line of the rectangle corresponds to the skin incision.

Skin incision: Under the microscope, approx. 25 mm long and approx. 40 mm from the midline. The surgeon's oculars are tilted 5°-10° towards the midline.

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Muscle splitting: With the index finger until contact is made with the junction of the medial third of the transverse process and the lateral joint facet. Alternatively: sequential dilators. The target area is exposed with narrow Langenbeck retractors.

Insertion of the spreading sleeve (SS): Advance a SS of suitable length between the Langenbecks with alternating rotating movements until contact is made with the corner between the lateral joint facet (medially) and the superior transverse process. Alternatively: after removing the Langenbecks, the SS is advanced through a dilator as far as the target area.

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Extraforaminal fenestration: Depending on the segment, particularly at L5/S1, partial burring of the lateral facet joint is indicated. Bipolar coagulation and incision of the medial insertion of the intertransverse muscle. Identification of the DH in the extraforaminal fat of the nerve. The latter is usually displaced by the DH in a cephalad and lateral direction as well as towards the surface. The extruded disc fragment is removed. If necessary, the IVS is partially cleared using standard microsurgical technique. The nerve root retractor can be replaced by the flattened sucker in right and left versions (3 sizes).

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The Dorsolateral Foraminotomy

Positioning: According to surgeon's preference: sitting, prone or so-called Concorde position (our preference).

X-ray labelling: (lateral view) Paramedian insertion of the cannula contralateral to the DH and projecting on the IVS. At the lower levels (C6/T1), the arms may be pulled with adhesive strapping or weights.

Skin incision: Under the microscope, approx. 18 mm long, 15 mm from the midline, centered on the lamino-facet junction. The surgeon's oculars are tilted 5°-10° towards the midline.

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Muscle splitting: The fascia and the single muscle layers are split with scissors or with a dissector according to their fibre direction. The lamino-facet junction is aimed at between two narrow Langenbecks. In this target area, the yellow ligament typically runs laterally in a V shape.

Insertion of the spreading sleeve (SS): Replace the Langenbecks with Love nerve root retractors and advance a SS of suitable length with alternating rotating movements until contact is made with bone. Alternatively: after removing the Langenbecks, advance the SS through a dilator as far as the target area. With both techniques, the SS is centered on the V-shaped yellow ligament.

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Interlaminar fenestration: Burr and/or punch out the lower rim of the superior hemilamina and/or the medial part of the joint complex (key-hole). Slim burr-handpieces facilitate visual control of the target structures. The yellow ligament is opened, the axilla of the nerve root is located, and the disc fragments are removed by standard microsurgical technique. The nerve root retractor can be replaced by the flattened sucker in right and left versions (3 sizes).

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